







## Eppendorf HeatSealer S100 and S200 Plate and Adapter Types

S100 (170 °C)	S200 (125 °C – 200 °C)	Suitable Plate Types
Plate Adapter, high profile 5391 070.034  	Plate Adapter, high profile 5392 070.011  	96-well Microplates 384-well Microplates  384-well PCR Plates  Deepwell Plates 384/200µL 96/500µL
96-well PCR Plate Adapter 5391 070.026  	96-well PCR Plate Adapter 5392 070.038  	96-well PCR Plates unskirted semi-skirted skirted
Plate Adapter, low profile 5391 070.018  	Plate Adapter, low profile 5392 070.020  	Deepwell Plates 96/1000µL 96/2000µL

Copyright© 2015 Eppendorf AG, Hamburg. All rights reserved, including graphics and images. No part of this publication may be reproduced without the prior permission of the copyright owner.

Eppendorf® and the Eppendorf logo are registered trademarks of Eppendorf AG, Germany.

Registered trademarks and protected trademarks are not marked in all cases with ® or ™ in this document.